











Совместимо с системой имплантатов

Совместимо с системой имплантатов










Astra Tech Osseospeed TX




*set (Scanbody/Multibase/CMFit Abutment + screw)






SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.27Hex 15Ncm ∅3.0 20Ncm ∅3.5/4.0 25Ncm ∅4.5/5.0

Astra Tech Osseospeed EV

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.27Hex 25Ncm

Description-Library code	1301- AS3	1302- ASR	1303- ASW
Connection Platform ∅	Yellow  ∅3.0	Aqua  ∅3.5/4.0	Lilac  ∅4.5/5.0
Product			
Scanbody	SC-AS3	SC-ASR	SC-ASW
GMS-Scanbody(Short)	-	GMS-ASRS	GMS-ASWS
GMS-Scanbody(Long)	-	GMS-ASRL	GMS-ASWL
LL-Multibase (Single)	GH: 1.3	LL2-AS3-T2	LL2-ASR-H2
	GH: 2.0	-	LL2-ASR20-H2
LL-Multibase (Bridge)	GH: 1.3	-	LL2-ASR-N2
	GH: 2.0	-	LL2-ASR20-N2
LS-Multibase / Single	LS2-AS3-T2	LS2-ASR-H2	LS2-ASW-H2
CMFit Abutment / ∅10	CMFit10-AS3-T2	CMFit10-ASR-H2	CMFit10-ASW-H2
CMFit Abutment / ∅14	-	CMFit14-ASR-H2	CMFit14-ASW-H2
3D-Analog	-	3D-ASR	3D-ASW
Screw	AS3C	ASRC	ASWC

Description-Library code	1304 - AEV30	1305 - AEV36	1306 - AEV42	1307 - AEV48	1308 - AEV54
Connection Platform ∅	Green  ∅3.0	Purple  ∅3.6	Yellow  ∅4.2	Blue  ∅4.8	Brown  ∅5.4
Product					
SC-Scanbody	SC-AEV30	SC-AEV36	SC-AEV42	SC-AEV48	SC-AEV54
GMS-Scanbody(Short)	-	-	-	-	-
GMS-Scanbody(Long)	-	-	-	-	-
LL-Multibase / Single	LL2-AEV30-T2	LL2-AEV36-T2	LL2-AEV42-T2	LL2-AEV48-T2	LL2-AEV54-T2
LL-Multibase / Bridge	LL2-AEV30-N2	LL2-AEV36-N2	LL2-AEV42-N2	LL2-AEV48-N2	LL2-AEV54-N2
LS-Multibase / Single	LS2-AEV30-T2	LS2-AEV36-T2	LS2-AEV42-T2	LS2-AEV48-T2	LS2-AEV54-T2
CMFit Abutment / ∅10	CMFit10-AEV30-T2	CMFit10-AEV36-T2	CMFit10-AEV42-T2	CMFit10-AEV48-T2	CMFit10-AEV54-T2
CMFit Abutment / ∅14	-	CMFit14-AEV36-T2	CMFit14-AEV42-T2	CMFit14-AEV48-T2	CMFit14-AEV54-T2
3D-Analog	3D-AEV30	3D-AEV36	3D-AEV42	3D-AEV48	3D-AEV54
Screw	AEV30C	AEV36C	AEV42C	AEV4854C	AEV4854C