


Совместимо с сирона Ti-base



Sirona Multibase



Совместим с inCoris ZI meso
Улучшенный дизайн для создания профиля прорезывания
Циркониевые блоки для Sirona multibase

Колпачок для сканирования Sirona Multibase

Description




	SHADE S2	SHADE S3
	S-MEZO-S2	S-MEZO-S3

Scan Cap

Колпачок для сканирования (короткий)	Колпачок для сканирования (длинный)
	
SCSN-CAP	SCSNL-CAP



Osstem Implant TS

*set (Sirona Multibase + screw)

	Platform Ø	TS Mini  Ø 3.0/3.5	TS Standard  Ø 4.0/4.5/5.0	Driver (Ncm)	1.2Hex
	Connection				
	Single [GH:1.2mm]	SRN2-SUM-H2	SRN2-SURO-H2	20Ncm	Ø 3.0/3.5
Single [GH:2.6mm]	-	SRN2-SURO26-H2	30Ncm	Ø 4.0/4.5/5.0	

Osstem Implant US / Biomet 3i Osseotite 4.1

*set (Sirona Multibase + screw)

	Platform Ø	US Regular  Ø4.1	Driver (Ncm)	1.2Hex
	Connection			
Single		SRN2-EXR-H2	30Ncm	



Совместимо с сирона Ti-base

Sirona Multibase

Description



Dentium Implantium / SuperLine

*set (Sirona Multibase + screw)

	Platform Ø	universal  Ø3.7/4.0/4.5	Driver (Ncm)	1.27Hex
	Connection			
	Single [GH:1.2mm]	SRN2-DER-H2	25~30Ncm	
Single [GH:2.6mm]	SRN2-DER26-H2			



Megagen AnyRidge

*set (Sirona Multibase + screw)

	Platform Ø	Universal  Ø3.5/4.0/4.5	Driver (Ncm)	1.2Hex
	Connection			
Single		SRN2-MARR-H2	35Ncm	



Megagen AnyOne

*set (Sirona Multibase + screw)

	Platform Ø	Universal  Ø 3.5/4.0/4.5/5.0	Driver (Ncm)	1.2Hex
	Connection			
	Single [GH:1.2mm]	SRN2-DERA-H2	25Ncm	
Single [GH:2.6mm]	SRN2-DERA26-H2			




Neo Biotech IS

*set (Sirona Multibase + screw)

	Platform Ø	 Ø3.5/4.0/4.5	Driver (Ncm)	1.2Hex
	Connection			
	Single [GH:1.2mm]	SRN2-DERN-H2	30Ncm	
Single [GH:2.6mm]	SRN2-DERN26-H2			

Straumann Bone Level

*set (Sirona Multibase + screw)

	Platform Ø	NC  Ø3.3	RC  Ø 4.1/4.8	Driver (Ncm)	Torx
	Connection				
	Single	SRN2-BLM-R2	SRN2-BLR-R2	35Ncm	